

Syringes | HPLC Autosampler

Manufactured since 2004, the comprehensive range of L-MARK® Gold and L-MARK Syringes for HPLC systems are designed and tested to meet critical instrument parameters.

Features and benefits

- Plunger guide maintains optimal plunger alignment for extended plunger lifetime.
- Syringes available with polyethylene, PTFE and black PTFE tipped plungers.
- Titanium plungers for longer lifetime.
- Needle setting in the syringe barrel is protected by a PTFE barrier for reduced exposure of the adhesive to samples.
- The L-MARK Gold Syringes include a chemically inert needle for extremely low carryover.

Product specifications

Available for CTC Analytics and Thermo Scientific HPLC autosamplers.

Specifications	
Fixed needle syringes	
Accuracy and reproducibility	±1% (dispensed volume)
Temperature range	5-70°C



For more information about this product visit www.trajanscimed.com or contact techsupport@trajanscimed.com